ENGLISH Instruction for the healthcare institution/provider to fill in a Patient Implant Card


RemeOs™ Screw Patient Implant Card for CE Market

The primary purpose of the RemeOs™ Screw implant card is to enhance patient safety by providing essential information about the implanted device, enabling quick identification in case of medical emergencies, recalls, or device-related issues.

A Patient Implant Card is to be filled in by the healthcare institution/provider. In the blank fields of the Patient Implant Card, fill in the following information:

1. Name of the patient or patient ID (marked with )
2. Date of implantation (marked with )
3. Name and address of the healthcare institution/provider which performed the implantation (marked with )

In the grey area of the Patient Implant Card, affix one sticker from the Patient Record Label that contains the device name, lot number, unique device identification in AIDC format, name, address, and the website of the manufacturer, use-by-date:

[Figure 1. Patient Record Label]

After filling in the information in the Patient Implant Card, give it to a patient.

Instruct the patient to keep the Patient Implant Card, for example, in his/her wallet. **This card must be carried at all times and presented to healthcare providers during medical consultations or emergencies** so that healthcare providers can identify the particular medical device implanted, follow the correct MR conditions, and use MRI safely.

MD

RemeOs™ Screw – Patient Implant Card

en

Device type: Absorbable bone screw

Affix Patient Record Label here

Bioretec Ltd. +358207789500

www.bioretec.com/remeos/mri-safety

Rev 1.0 2025-01-13, 109129

Front

Non-clinical testing demonstrated that the RemeOs™ Screw is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury:	
Static Magnetic Field Strength (B ₀)	1.5 T or 3.0 T
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg
Scan Duration	Under the scan conditions defined, the implant can be scanned continuously for 25 minutes. With 5 minutes break, this can be repeated twice in 60 minutes.

Back

The outer dimensions are the same as those of the credit cards: Dimensions: 85.60 x 53.98 mm, radius 2.88-3.48 mm³. Text at least 2 mm high.

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Bioretec Ltd. +358207789500 www.bioretec.com/remeos/mri-safety