ENGLISH

PATIENT INSTRUCTIONS: POSTOPERATIVE CARE

Appropriate immobilization must be used, and the patient must be detailed instructed for postoperative care because premature bending, loosening, fracture or migration of the implant may result from too early weight bearing or physical activity.

Absorption and Healing Period

You have undergone an operation where an absorbable bone fixation device (RemeOs™ Screw PG) manufactured by Bioretec Ltd. has been implanted in your body.

This bone fixation device is made of material that corrodes gradually due to interaction with body fluids into magnesium, calcium and zinc oxides and hydroxides, which are naturally occurring elements in human body and are known to promote bone regeneration.

Foreign materials (e.g., non-absorbable implants) could cause harm over the years. However, after absorption of the RemeOs™ Screw PG, the tissue is free from foreign material. Thus, there is no need for removal surgery for a bone fixation device after the bone tissue healing.

The strength of the device is maintained for at least 8 weeks, which is an adequate period for bone to heal. The final absorption takes place within approximately 2-3 years after the surgery.

Magnetic Imaging Conditions

After the fixation device is absorbed, it does not disturb magnetic imaging. During the absorption period, the fixation device is magnetic imaging conditional. This means that during the absorption period, you can be scanned safely in an MR system under the following conditions:

Static Magnetic Field Strength (B0)	1.5T or 3.0T
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2W/kg
Scan Duration	Under the scan conditions defined, the implant can be scanned continuously for 25 minutes. With 5 minutes break, this can be repeated twice in 60 minutes.

Patient Implant Card

Keep your Patient Implant Card carefully with you and show the card to the staff when dealing with healthcare so that the healthcare staff can treat you correctly. The Patient Implant Card guides the healthcare personnel to use the correct conditions during the MR imaging.

MD

RemeOs™ Screw – Patient Implant Card

en

?

Device type: Absorbable bone screw

Affix Patient Record Label here

Bioretec Ltd. +358207789500 www.bioretec.com/remeos/mri-safety

Rev. 1.0 2025-01-13, 109129

Front

Non-clinical testing demonstrated that the RemeOs™ Screw is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions. Failure to follow these conditions may result in injury:

Static Magnetic Field Strength (B ₀)	1.5 T or 3.0 T
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg
Scan Duration	Under the scan conditions defined, the implant can be scanned continuously for 25 minutes. With 5 minutes break, this can be repeated twice in 60 minutes.

Back

Adverse Effects

Reported adverse effects for absorbable bone fixation devices are tissue reactions towards implant material. These show typically as redness, swelling, and fluid accumulation in the operated area. Contact your doctor if you have any concerns about your adverse effects.

Postoperative Reminders

The following information describes basic after-care recommendations that apply to many types of surgery. These basic suggestions to improve your recovery may be supplemented by other directions from your doctor and nurse.

For Same Day Surgery patients, please refer to the following guidelines:

Activity:

- Make arrangements to have a ride home after your surgery and to have someone stay with you when you arrive home as long as your doctor has instructed.
- Plan to rest and relax at least until the evening of the next day following your operation. Your physician will give you more detailed instructions.
- Check with your physician before performing any strenuous activities after the surgery.

Medications

- Your doctor will most likely prescribe medications during your recovery period. Please have the prescriptions filled as soon as you leave the hospital.
- Take the medicine as instructed for the maximum benefit.
- Follow doctor`s detailed instructions for postoperative care, including weight bearing limitations and physical activity.

In case of an increase in redness, pain, swelling, or fluid accumulation in the operated area, you must get in contact with your doctor immediately to find out the cause of the illness and for the treatment.

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